

MS08 H₂O₂ Amperometry

Portable and laboratory measuring instrument



Multiparameter system

Check for the presence of hydrogen peroxide

Hydrogen peroxide finds many applications in various fields such as medicine and the food industry which use it for disinfection and sterilization, but also in the manufacture of cosmetics, the paper industry and in particular in the treatment of the water. H₂O₂ removes organic pollutants. However, it is also involved in corrosion phenomena and the generation of oxidizing species that are very aggressive for the body. Controlling the dissolved H₂O₂ concentration is an important parameter for adapting water treatment, especially for the oxidation of H₂S.

The precise and reliable in situ determination of the concentrations of this dissolved gas is possible with the MS08-H₂O₂ for online or portable measurements. The micro H₂O₂ sensor allows rapid measurement with very high resolution. The multi-parameter MS08 system collects raw H₂O₂ information and performs compensation calculations based on the temperature of the medium.

Applications

- Monitoring and protection of wastewater networks
- Control of H₂O₂ injections
- Industrial process management
- Protection of personnel before intervention

Advantages

- Measurement without direct sampling in the medium
- No interference with turbidity or water color
- Automatic temperature compensation
- Versatile - portable / laboratory and stationary measurements
- Direct display of measurement data
- Data extraction to computer by USB

Online / portable measurement

The MS08 box is designed for portable measurements (approx. 12 hrs battery life) and can also be connected to a 220VAC mains socket for continuous measurements.




The dissolved H₂O₂ passes through the gas permeable silicone membrane. It diffuses to the working electrode where an electrochemical oxidation reaction takes place. The current generated, proportional to the concentration of hydrogen peroxide, is measured by the probe.

This current from 0 to 400 pico-amperes is then converted and exploited by the MS08 box, the measurement data are compensated using the temperature measurement.

MS08 H₂O₂ Amperometry

Technical specifications

Measuring principle	Amperometric measurement
Technology	Micro membrane sensor with redox catalyst
Temperature compensation	Automatic Pt100, Pt1000
Electrical polarization time	Automatic from 5 to 15 minutes maximum wait at start-up
Measuring ranges	Type I 0,02...10 %H ₂ O ₂ Type S Other ranges on request
Response time	T 90% < 2 secondes
Measuring principle	1% of the measured value
H ₂ O ₂ consumption	Negligible
Working pH range	from 0 to 11 pH
Probe body materials	H ₂ O ₂ probe - titanium / temperature probe - plastic
Dimensions (d x L)	H ₂ O ₂ probe - 17 mm x 205 mm / temperature-pH probe - 12 mm x 120 mm
MS08 power supply	6 Mignon type batteries / 220 VAC with charger supplied
Interface	Digital Data display on the integrated screen - H ₂ O ₂ / T ° / residual current / pH (option) RS-232 / USB (option)
	Analog 4 ... 20 mA on request
Data acquisition	On computer, software not supplied, downloadable free of charge type «HYPERTERMINAL»
Transmission / frequency	ASCII string / 2 seconds
Operation	Conversion by software for access to Lotus 1-2-3 or Excel type calculation software
Service life of the micro H ₂ O ₂ sensor	6 months in portable use, 10 months continuously (depends on stress due to variations in H ₂ O ₂)
Interference on measurement	No interference in salt water up to 40 g / l salt No interference in the presence of: carbon dioxide, oxygen, methane, hydrogen, ammonia, carbon monoxide, organic solvents (maximum 20% vol.), Acetic acid, sulphide dimethyl, HCN and solid compounds
Maintenance	Cleaning the measuring membrane with distilled water after each use
Temperature of the medium / sample	0 ... + 30 ° C (40 ° C possible with a specific calibration on request)
Ambient temperature	0 ... + 40 ° C
Storage temperature	0 ... + 40 ° C



Combined pH-T